


# Spatially Resolved Layer Thickness Analysis of Thin Alloy Coatings with Micro-XRF



Bruker Nano Analytics, Berlin, Germany  
Webinar, May 12, 2016



# Presenters



Falk Reinhardt

Application Scientists micro-XRF,  
Bruker Nano Analytics, Berlin, Germany



Samuel Scheller

Sr. Product Manager micro-XRF,  
Bruker Nano Analytics, Berlin, Germany

# Overview



- Introduction
- The M4 Tornado micro-XRF Instrument
- Example 1 – Ti / ZnAlMg coating on sheet metal
- Example 2 – Cu-Al coating on glass
- Summary
- Live demonstration
- Questions

# Introduction



- Please also have a look at our Webinar on Layer Analysis:  
"Automated Analysis of Decorative and Functional Coatings with the M1 MISTRAL Micro-XRF Instrument" from the 4<sup>th</sup> May 2016
- The M4 Tornado is a high speed micro-XRF instrument for challenging applications and element distribution analysis
- We want to present you with a couple of analysis examples looking at metal and alloy layers,
- Also, we want to provide a live demonstration of the instrument
- Any questions you may have, we want to answer and discuss

# M4 Tornado Micro-XRF Spectrometer

## Standard configuration



### **30 W microfocus X-ray Tube**

- Rh target
- Polycapillary lens with spot size  $< 20 \mu\text{m}$
- Resolve micro structures
- Effective for all material types

### **Silicon drift detectors (SDD)**

- $30 \text{ mm}^2$  with energy resolution  $< 145 \text{ eV}$  (for Mn-K $\alpha$ )
- Optional  $60 \text{ mm}^2$  with energy resolution  $< 145 \text{ eV}$
- Second detector option for double speed

# M4 Tornado Micro-XRF Spectrometer

## Standard configuration



### **Vacuum sample chamber**


- Pressure range from 1 mbar up to atmosphere
- Light element detection down to Sodium
- Stage 190 mm x 160 mm
- Maximum sample height of 120 mm
- 5 kg capacity
- Stage measurement speed up to 100 mm/s
- 4  $\mu\text{m}$  minimum step size

### **Advantages**

- Non-destructive, spatially resolved material analysis
- Requires no sample preparation
- Quality control screening measurement of structured samples
- X-ray „depth view“ is ideal for analysis of thin alloy layers




# Position Tagged Spectrometry Imaging The HyperMap Database




# X-Ray Fluorescence ...in Coating Analysis



## Excitation



## XRF Detection




## Silicon **D**rift **D**etector with XFlash® Technology

- X-rays can penetrate and excite matter
- Signal from base material and covered layers can still be detected
- X-rays are attenuated in characteristic ways on their path through matter
- Intensity ratios of observed elemental lines, allows for calculation of respective layer thicknesses





# Example 1



Ti /ZnAlMg coating on sheet metal

# Ti /ZnAlMg: Coating on Sheet Metal

## Overview and Measurement Data (Map)



### Mapping parameter

Width:	1640	pixel
	164	mm
Height:	694	pixel
	69.4	mm
Pixel Size:	100	µm
Total number of pixel:	1138160	pixel

### Acquisition parameter



Frame count:	1
Pixel time:	15 ms/pixel
Measure time:	17072 s
Overall time:	17705 s

### Tube parameter


High voltage:	50	kV
Anode current:	600	µA
Filter:	Al100	
Optic:	Lens	
Chamber at:	Air 1023.3	mbar

# Ti /ZnAlMg: Coating on Sheet Metal

## Element Distribution



Measurements show small variations of the Zn and Fe intensities. Ti peak is relatively small, but looks to be distributed homogeneously in the sample.



## Example 2



Cu-Al coating on glass



# Cu-Al Coating on Glass

## Overview and Measurement Data (Map)




- Glass substrate is coated through Magnetron sputtering using a dual Cu-Al-target
- Layer thickness and gradient of Cu:Al ratio were to be measured
- Over night measurement (15 h) over the whole sample was performed
- 50  $\mu\text{m}$  step size, 50 ms dwell per pixel

Cu:Al gradient layer, manufactured through Magnetron sputtering of a dual Cu/Al target by K. Harbauer, W. Kohrt, K. Ellmer, Helmholtz-Zentrum Berlin für Materialien und Energie

# Cu-Al Coating on Glass

## Overview and Measurement Data (Map)






Map information			
<b>Mapping parameter</b>			
Width:	1000	pixel	
	50	mm	
Height:	1000	pixel	
	50	mm	
Pixel Size:	50	µm	
Total number of pixel:	1000000	pixel	
<b>Acquisition parameter</b>			
Frame count:	1		
Pixel time:	50	ms/pixel	
Measure time:	13:53 h		
Overall time:	14:58 h		
<b>Tube parameter</b>			
High voltage:	50	kV	
Anode current:	200	µA	
Filter:	Empty		
Optic:	Lens		
Chamber at:	Vacuum 20	mbar	
Anode:	Rh		
<b>Detector parameter</b>			
Selected detectors:	1		
Close			

Cu:Al gradient layer, manufactured through Magnetron sputtering of a dual Cu/Al target by K. Harbauer, W. Kohrt, K. Ellmer, Helmholtz-Zentrum Berlin für Materialien und Energie



# Cu-Al Coating on Glass

## Element Distribution




Qualitative distribution of the major elements comprised of base and layer

- Copper gradient distribution from top-right corner to bottom left
- The maximum of the Aluminum distribution seems offset to the left corner

# Cu-Al Coating on Glass

## Layer Thickness Analysis



No standards available

Hence pure FP Quantification

5x5-Binning:

- effectively 1.25 s per Pixel
- 250  $\mu\text{m}$  spatial resolution


The layer is not homogeneous. The center of the radially symmetric thickness gradient is slightly eccentric. The thickness of the Cu:Al layer reduces from  $\sim 380$  nm at the center down to  $\sim 250$  nm at the edges.

Structure		
Layer	Chemical elements	
layer1	Cu	Al
base	Si	

Measure conditions	
Parameter	Value
HV / kV	50
Current / $\mu\text{A}$	200
Collimator	25 $\mu\text{m}$ LENS
Atmosphere	Vacuum

# Cu-Al Coating on Glass

## Comparison (Point Measurements)




Alternatively, to reduce the measurement time, an array of point measurements can be performed. The inserting of a grid of measurement point is an automatic function in the software under Multi-Point. Left shows a 10x10 grid. Measurement time is 10 s per point.

Each spectrum can be quantified by the same layer analysis method. The results (layer thicknesses and composition) are shown in tabular form and can be exported to Excel in the .xls format. This allows, for example, to visualize the result of the element and layer thickness distribution.

# Cu-Al Coating on Glass

## Layer Thicknesses and Composition



Layer thickness

0.239	0.259	0.291	0.301	0.311	0.321	0.325	0.314	0.000	0.000
0.258	0.287	0.307	0.318	0.346	0.340	0.350	0.000	0.000	0.300
0.269	0.298	0.322	0.342	0.355	0.361	0.000	0.000	0.332	0.316
0.277	0.314	0.324	0.347	0.369	0.366	0.000	0.367	0.347	0.327
0.284	0.325	0.342	0.359	0.378	0.000	0.372	0.370	0.355	0.331
0.283	0.312	0.346	0.360	0.000	0.379	0.386	0.367	0.357	0.324
0.281	0.314	0.338	0.000	0.301	0.382	0.369	0.370	0.355	0.309
0.264	0.294	0.000	0.000	0.350	0.356	0.370	0.347	0.337	0.305
0.249	0.000	0.030	0.319	0.326	0.342	0.336	0.327	0.319	0.291
0.000	0.000	0.274	0.289	0.295	0.303	0.302	0.295	0.286	0.256

Cu Concentration

77.2	79.0	79.3	81.3	82.4	83.2	84.5	85.6	0.0	0.0
76.6	77.4	78.6	80.4	80.1	82.3	82.6	0.0	0.0	86.3
75.4	76.9	77.7	78.6	80.0	80.8	0.0	0.0	85.2	85.1
74.4	74.4	77.5	78.1	78.3	80.0	0.0	82.2	83.6	84.0
72.9	73.0	74.7	75.8	76.4	0.0	79.6	80.4	81.6	83.0
71.3	72.6	72.7	74.5	0.0	76.4	77.0	79.5	80.5	82.2
69.7	70.4	71.2	0.0	74.7	74.5	76.9	77.3	79.3	81.8
69.3	69.7	0.0	0.0	72.9	74.3	74.2	77.1	78.2	79.8
67.5	0.0	0.0	70.4	71.9	72.8	74.6	76.3	77.2	78.4
0.0	0.0	68.2	69.4	71.4	72.7	74.3	75.8	76.9	79.0

Al Concentration

22.8	21.0	20.7	18.7	17.6	16.8	15.5	14.4	0.0	0.0
23.4	22.6	21.4	19.6	19.9	17.7	17.4	0.0	0.0	13.7
24.6	23.1	22.3	21.4	20.0	19.2	0.0	0.0	14.8	14.9
25.6	25.6	22.5	21.9	21.7	20.0	0.0	17.8	16.4	16.0
27.1	27.0	25.3	24.2	23.6	0.0	20.4	19.6	18.4	17.0
28.7	27.4	27.3	25.5	0.0	23.6	23.0	20.5	19.5	17.8
30.3	29.6	28.8	0.0	25.3	25.5	23.1	22.7	20.7	18.2
30.7	30.3	0.0	0.0	27.1	25.7	25.8	22.9	21.8	20.2
32.5	0.0	0.0	29.6	28.1	27.2	25.4	23.7	22.8	21.6
0.0	0.0	31.8	30.6	28.6	27.3	25.7	24.2	23.1	21.0



# Cu-Al Coating on Glass

## Results Comparison




22.8	21.0	20.7	18.7	17.6	16.8	15.5	14.4	0.0	0.0
23.4	22.6	21.4	19.6	19.9	17.7	17.4	0.0	0.0	13.7
24.6	23.1	22.3	21.4	20.0	19.2	0.0	0.0	14.8	14.9
25.6	25.6	22.5	21.9	21.7	20.0	0.0	17.8	16.4	16.0
27.1	27.0	25.3	24.2	23.6	0.0	20.4	19.6	18.4	17.0
28.7	27.4	27.3	25.5	0.0	23.6	23.0	20.5	19.5	17.8
30.3	29.6	28.8	0.0	25.3	25.5	23.1	22.7	20.7	18.2
30.7	30.3	0.0	0.0	27.1	25.7	25.8	22.9	21.8	20.2
32.5	0.0	0.0	29.6	28.1	27.2	25.4	23.7	22.8	21.6
0.0	0.0	31.8	30.6	28.6	27.3	25.7	24.2	23.1	21.0

**Element distribution of Al** (qualitative, left) shows maximum is not in the corner, the layer composition (quantitative, right) shows that the highest Al:Cu ratio is in the corner. The layer thicknesses is also thinner there.



77.2	79.0	79.3	81.3	82.4	83.2	84.5	85.6	0.0	0.0
76.6	77.4	78.6	80.4	80.1	82.3	82.6	0.0	0.0	86.3
75.4	76.9	77.7	78.6	80.0	80.8	0.0	0.0	85.2	85.1
74.4	74.4	77.5	78.1	78.3	80.0	0.0	82.2	83.6	84.0
72.9	73.0	74.7	75.8	76.4	0.0	79.6	80.4	81.6	83.0
71.3	72.6	72.7	74.5	0.0	76.4	77.0	79.5	80.5	82.2
69.7	70.4	71.2	0.0	74.7	74.5	76.9	77.3	79.3	81.8
69.3	69.7	0.0	0.0	72.9	74.3	74.2	77.1	78.2	79.8
67.5	0.0	0.0	70.4	71.9	72.8	74.6	76.3	77.2	78.4
0.0	0.0	68.2	69.4	71.4	72.7	74.3	75.8	76.9	79.0

**Element distribution of Cu** (qualitative, left) The layer composition (quantitative, right) shows the highest Cu:Al ratio in the corner.



0.299	0.299	0.291	0.301	0.311	0.321	0.325	0.314	0.000	0.000
0.258	0.287	0.307	0.318	0.346	0.340	0.350	0.000	0.000	0.300
0.269	0.298	0.322	0.342	0.355	0.361	0.000	0.000	0.332	0.316
0.277	0.314	0.324	0.347	0.369	0.366	0.000	0.367	0.347	0.327
0.284	0.325	0.342	0.359	0.378	0.000	0.372	0.370	0.355	0.331
0.283	0.312	0.346	0.360	0.000	0.379	0.386	0.367	0.357	0.324
0.281	0.314	0.338	0.000	0.301	0.382	0.369	0.370	0.355	0.309
0.264	0.294	0.000	0.000	0.350	0.356	0.370	0.347	0.337	0.305
0.249	0.000	0.030	0.319	0.326	0.342	0.336	0.327	0.319	0.291
0.000	0.000	0.274	0.289	0.296	0.303	0.302	0.295	0.286	0.256

Very similar results with respect to position and min/max-values of the **layer thickness distribution**.

Spatial resolution vs. measurement time (here 250  $\mu\text{m}$  with 15 h vs. 5 mm with < 20 min)

# Summary



- Micro-XRF allows measurement of thin (micrometer→nanometer) layer thicknesses with high ( $< 20 \mu\text{m}$ ) lateral spatial resolution
- Through the HyperMap database, element distributions can be saved with spectra and processed afterwards
- Intensity variations are not necessarily proportional to the layer thickness changes
- To determine the thickness of a layer, it must be quantified or calculated by applying a layer thickness method
- This can then be visualized as distribution or layer thickness image
- Alternatively, using single point measurements, the layer can be analyzed and quantified, followed by additional processing in Excel to provide a visual representation of the element and layer thickness distribution



## Any Questions?

Please type in the questions you might have  
in the Q&A box and press Submit.


**For more information, please contact us:**

[falk.reinhardt@bruker.com](mailto:falk.reinhardt@bruker.com)

[samuel.scheller@bruker.com](mailto:samuel.scheller@bruker.com)

[info.bna@bruker.com](mailto:info.bna@bruker.com)

<https://www.bruker.com/m4tornado>



Innovation with Integrity